



 **LASANI
SURGICAL**

EDITION 2026

PODO
COLLECTION DE BEAUTY



Face



Hand



Foot



Designed by
LASANI Surgical Industry
PAKISTAN




TRADITION & INNOVATIVE TECHNIQUE





PODLOGIE

& BEAUTY CARE INSTRUMENTS





Welcome to LASANI SURGICAL Instruments

Since our inception, our unwavering dedication has centered around providing professionals worldwide with premium instruments at optimal value. Specializing in the distribution of Orthodontics, General Dental, Podiatry, and Surgical instruments

LASANI SURGICAL distinguishes itself through superior quality and a commitment to surpassing industry benchmarks. Our instruments embody excellence, crafted through a synthesis of advanced technology, precision-driven CNC machines, and the expertise of skilled mechanical engineers. From the initial point of contact through customer service, our focus spans the entire spectrum of operations, consistently aiming to build enduring relationships and exceed expectations.

At LASANI SURGICAL, you are not merely a customer; you are our esteemed partner, benefiting from our unwavering commitment to professionalism and elegance in every facet of our service. We are grateful and thank you for partnering with LASANI SURGICAL. We truly value you and your business and look

Sincerely,

KABEER HUSAIN

LASANI SURGICAL INSTRUMENTS





OUR MISSION



Mission

Lasani Surgical Instruments Mission

It is our mission to lead the global industry by providing the highest quality instruments at the best value through dependable, personalized service, on-time delivery, and continuous improvement of our instruments and processes. LASANI SURGICAL stands behind every instrument with a lifetime guarantee.



Vision

Lasani Surgical Instruments Vision

Our vision is to be a leading medical instruments supplying company in the world, expanding our services to reach additional peoples. We work to be trusted by customers.



Values



INNOVATION



CLASS



QUALITY

BEAUTY
DÜSSELDORF

31.3.–2.4.2023
www.beauty.de



"Stepping into Success at Podo Instruments:

Where Every Stride Marks a Triumph!

Amidst the symphony of innovation and expertise, our presence at Podo Instruments echoed triumphantly. The stage was set, the spotlight shone on our breakthrough products, and the applause reverberated through the halls, a testament to our resounding success. As we showcased our cutting-edge offerings and engaged in insightful conversations, our journey through Podo Instruments became a defining chapter in our pursuit of excellence."





PACKAGING



Trimming



Forgings



Milling



Fitting & Grinding



Finishing

Lasani Surgical Instruments **Manufacturing Process**

Here we will shortly brief you about the process of manufacturing



Tradition in Line with Innovation

Small masterpieces made from steel, which also meet the highest professional standards and prove themselves through their function, aesthetics, ergonomics and long-term precision cutting ability.

In addition, the product properties such as hardness, robustness and corrosion resistance significantly determine the performance, reliability and durability of the instruments. From the very first look, excellent craftsmanship and eye catching appearance of our instruments produces the right impression. Using of only the best German stainless steel and precise finishing processes insures that our instruments are made to last and highly appreciated by beauty professionals all over the world





WE BRING THEM TO LIFE










TABLE OF CONTENTS

16

NAIL & INGROWN NIPPERS

/ NAGELZANGEN & ZANGEN FÜR EINGEWACHSENE NÄGEL

- Head Cutters 14cm (Page no 18)
- Pedicure Nail Nippers 14cm (Page no 19-21)
- Head Cutter 12cm (Page no 22)
- Nail Nipper Curved Lap Joint D/S 13cm (Page no 22)
- Pedicure Nail Nippers 12cm (Page no 23)
- Double Action Nail Nippers 15cm (Page no 24)
- Double Action Nail Nippers 12cm (Page no 25)
- Double Action Nail Nipper Straight Fine 15cm (Page no 25)
- Nail Splitter 13cm (Page no 26)
- Nail Pulling Forcep 14cm (Page no 26)
- Nail Nipper Curved Lap Joint S/S 12cm & 10.5cm (Page no 27)
- Toenail Nipper Curved Jaw Barrel Spring 12cm & 10.5cm (Page no 27)
- Ingrown Nail Nippers 10.5cm (Page no 29)
- Ingrown Nail Nippers 11.5cm (Page no 30-31)
- Ingrown Nail Nippers 13cm (Page no 32-33)
- Ingrown Nail Nippers 15cm (Page no 34)
- Podology Set (Page no 35)
- Ingrown Nail Nippers Scissor Ring (Page no 36)
- Nail Clippers (Page no 37)
- Forming and Cutting Pliers (Page no 38-40)



OUR VISION

Our vision is to be one of the leading instrument company in the world, expanding our services to reach additional peoples. We work to be trusted by customers.

41 CUTICLE NIPPERS / HAUTZANGEN

- Cuticle Nippers 10cm D/S (page no 42)
- Cuticle Nippers 10cm S/S (page no 43)
- Cuticle Nippers 11.5cm D/S (page no 44)
- Cuticle Nippers 11.5cm S/S (page no 45)
- Cuticle Nippers 13cm D/S (page no 46)
- Cuticle Nippers 13cm S/S (page no 47)
- Cuticle Nippers 11.5cm D/S with lock (page no 48)
- Acrylic Nail Nipper 11.5cm S/S with lock (page no 49)
- TIP Slicer with Acrylic nail set (page no 50)
- Cuticle Nippers Scissor Ring (page no 51)
- Cuticle Nippers 12cm & 13cm lock (page no 51)
- Cuticle Nippers 10cm Lap Joint (page no 52)

53 CUTICLE PUSHERS / NAGELHAUTDRÜCKER

- Cuticle Pushers Sharp (page no 54-56)
- Cuticle Flat Pushers (page no 54-56)
- Cuticle Remover Pushers (page no 54-56)
- Cuticle Round Pushers (page no 54-56)
- Cuticle Pushers with Knife (page no 54-56)
- And More.... (page no 54-56)





57 EXCAVATORS, CURETTES & PROBES

- Excavators double sided (page no 58-60)
- Curettes double sided (page no 58-60)
- Probes Single sided (page no 58-60)
- Chisels Single sided (page no 58-60)

61 EXTRACTORS

- Extractors double sided (page no 62-63)
- Extractors Single sided (page no 62-63)
- Extractors Tweezer (page no 63)
- Chisel Handles (page no 64)
- Scalpel Handles (page no 64)





65 INGROWN NAIL FILES / NAGELFEILEN FÜR EINGEWACHSENE NÄGEL

- Ingrown Nail Files Double Sided (page no 66-71)
- Ingrown Nail Files Single Sided (page no 66-71)

72 SPATULAS & NAIL PACKERS

- Spatulas Double Sided (page no 73-75)
- Spatulas Single Sided (page no 73-75)
- Nail Packers (page no 73-74)
- Mirror (page no 75)

76 EYEBROW & EYELASH TWEEZERS / PINZETTEN

- Eyebrow Tweezers Standard Tip (page no 77-78)
- Eyebrow Tweezers Pattern/colourfull (page no 79)
- Eyebrow Tweezers Narrow Tip (page no 80)
- Body Waxing Tweezer (page no 80)
- Eyelash Tweezers (page no 82-84)
- Tissue Forceps Tweezers (page no 85)

86 SCISSORS / SCHEREN



- Stork scissors (page no 87-92)
- Curved nail scissors (page no 87-92)
- Toenail scissors (page no 87-92)
- Fiberglass scissors (page no 87-92)
- Ultra fine scissors (page no 87-92)
- Mustache scissors (page no 87-92)
- Safety scissors (page no 87-92)
- Cuticle scissors fine/curved,pointed (page no 87-92)
- And More..... (page no 87-92)
- Needle Holders (page no 94)
- Mosquito Forceps (page no 95)
- Mathieu Needle Holder (page no 95)





96 FOOT FILES/RASPS, NAIL FILES ,TRAYS & KITS

- Stainless steel foot files (page no 96-97)
- Wooden foot files (page no 98-99)
- Plastic foot files (page no 99)
- Foot rasp & files (page no 100-101)
- Stainless Steel nail files (page no 102-103)
- Instruments trays (page no 104-105)
- Instruments kits (page no 106-107)



www.lasanisurgicalind.com

CONTACT US

Gerhard Stalling Str. 70,
26135 Oldenburg, Germany
+49 441 40579913

info@lemedinst.com

info@lemedinstruments.de



NAIL NIPPERS

LASANI SURGICAL Nail nippers are manufactured using high-pressure forging of layers made of 100% surgical-grade stainless steel, capable of performing under normal or expected conditions and achieving Rockwell HRC 48–50. Our nippers are designed for professionals for durability, sharp cutting, and long-lasting

The nail nippers are formed with a box-joint method which prevents lateral play or staggering and retains cutting-edge alignment for a lifetime. Both blades are

Our nail nippers' curved jaws are designed to allow a clean cut and provide a slightly curved nail shape without ripping or tearing. Blades touch evenly when closed, and both blades are exactly aligned to cut hard nails efficiently. Double springs are inserted to allow smooth and precise floating movement

The handles have a round tension to give you the highest possible grip which allows the hand to have maximum control.


Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, it features gently curved, palm-friendly, easy-to-grip ergonomically designed grooved finished handles which gives optimal comfort, performance and expresses the ultimate in functional beauty.

LASANI SURGICAL Nail nippers come in various forms such as nail cutters, head cutters, and ingrown nail cutters. Jaw styles include curved, straight, flat, ultra-fine, and pointed to maximize comfort and user-


All LASANI SURGICAL nippers have springs assembled with a screw that can



Available In Two Handles
Groove Handle



Plain Handle



High Quality
Instruments

INNOVATION AND RELIABLE QUALITY

Front Working Jaw

Nail Nipper's jaw are designed for perfect work.

Sharpness

Long Lasting Sharpness. Box-joint Nail Nipper have Quality to stops side-ways movement and It adjusts the alignment of the cutting edges.

Hand Crafted Box Joint Construction

Fully Tempered Box locks maintain instruments alignment. The most reliable design possible

Lifetime Warranty

Against forging breakage

Stainless Steel

Forging is made with highest quality durable (**Surgical Grade Stainless Steel**).

Design & Performance

Outstanding instrument's Design, size and functions for optimal comfort and performance.

Springs assembled with a screw which can be replaced.
Satin Finish simplifies easy cleaning.

Drawings



Lap-joint



Box-joint





LSI-00-201

Head Cutter

Double Spring with Lock

Size: 14cm

Jaw Size: 11mm



- Fully Tempered Box locks maintain instruments alignment. Box-joint Nail Nipper have Quality to stops side-ways movement and It adjusts the alignment of the cutting edges.
- Nail Nipper's jaw are designed for perfect work. Every sharp cutting. Long Lasting Sharpness. The quality of the steel has a crucial influence on the quality of Product. It is fundamental for the Sharpness and durability of the product
- Springs assembled with a screw which can be replaced. Outstanding instrument's Design, size and functions for optimal comfort and performance.
- Lock work Smoothly. Satin Finish simplifies easy cleaning.




LSI-00-202

Head Cutter

Double Spring with Lock

Size: 14cm

Jaw Size: 22-23mm




LSI-00-203

Head Cutter

Double Spring with Lock

Size: 14cm

Jaw Size: 20mm





LSI-00-204
Pedicure Nail Nipper
Straight
Double Spring with Lock
Size: 14cm




LSI-00-205
Pedicure Nail Nipper
Curved
Double Spring with Lock
Size: 14cm






LSI-00-206
Pedicure Nail Nipper
Angular
Double Spring with Lock
Size: 13.5cm




LSI-00-207
Pedicure Nail Nipper
Extra-Fine Pointed
Double Spring with Lock
Size: 13.5cm






LSI-00-208
Pedicure Ingrown Nail Nipper
Double Spring with Lock
Size: 13.5cm



LSI-00-209
Pedicure Nail Nipper Fine Pointed Jaw
Double Spring with Lock
Size: 13.5cm



Pedicure Nail Nippers



Box-joint



HEAD CUTTER

5.5" inch jaw



LSI-00-210

Head Cutter

Double Spring with Lock

Size: 12cm

Jaw Size: 17mm



Box-joint

PREMIUM QUALITY

LSI-00-211


Nail Nipper Curved

Double Spring

Size: 13cm




Lap-joint




LSI-00-212
Pedicure Nail Nipper
Straight
Double Spring with Lock
Size: 12cm



LSI-00-213
Pedicure Nail Nipper
Curved
Double Spring with Lock
Size: 12cm



LSI-00-214
Pedicure Ingrown
Nail Nipper
Double Spring with Lock
Size: 12cm




LSI-00-215
Pedicure Nail Nipper
Extra-Fine Pointed
Double Spring with Lock
Size: 12cm



Pedicure Nail Nippers





Box-joint




LSI-00-216
Double Action
Nail Nippers Curved
Wire Spring with Lock
Size: 15cm




LSI-00-217
Double Action
Nail Nippers Straight
Wire Spring with Lock
Size: 15cm








LSI-00-218
Double Action
Nail Nippers Curved
Wire Spring with Lock
Size: 12cm



LSI-00-219
Double Action
Nail Nippers Straight
Wire Spring with Lock
Size: 12cm




LSI-00-220
Nail Nippers Double
Action Curved Fine
Double Spring
Size: 15cm




LSI-00-226
Nail Nippers Double
Action Straight Fine
Double Spring
Size: 15cm




Double Action Nail Nippers



LSI-00-2227
Platypus Nail Pulling
Forceps Narrow
5-1/2"
Size: 14cm




LSI-00-2223
Nail Splitter
Standard Jaw
Double Spring
Size: 13cm





LSI-00-224
Nail Nipper Curved
Single Spring
Size: 12cm



LSI-00-225
Toenail Nipper
Curved Jaw
Barrel Spring with Lock
Size: 13cm



LSI-00-226
Nail Nipper Curved
Single Spring
Size: 10.5cm



LSI-00-227
Toenail Nipper
Curved Jaw
Barrel Spring with Lock
Size: 12cm






INNOVATION WITH PERFECTION



INGROWN NAIL NIPPERS

LASANI SURGICAL Ingrown Nail Nippers come in 4 unique handle sizes, making the nippers user-friendly. The nail nippers are formed using a box-joint method that prevents lateral play or staggering and retains cutting-edge alignment for a lifetime. Both blades are precisely aligned in proper relative positions to cut hard and ingrown nails efficiently.



LSI-00-228
Ingrown Nail Nipper
Pointed Jaw
Double Spring
Size: 10.5cm



LSI-00-229
Ingrown Nail Nipper
Ultra Fine
Pointed Jaw
Double Spring
Size: 10.5cm



LSI-00-230
Ingrown Nail Nipper
Standard Jaw
Double Spring
Size: 10.5cm



Ingrown Nail Nippers



LSI-00-231
Ingrown Nail Nipper
Standard Jaw
Double Spring
Size: 11.5cm



LSI-00-232
Ingrown Nail Nipper
Pointed Jaw
Double Spring
Size: 11.5cm



LSI-00-233
Ingrown Nail Nipper
Slim Jaw
Double Spring
Size: 11.5cm



LSI-00-234
Ingrown Nail Splitter
Double Spring
Size: 11.5cm






LSI-00-235
Ingrown Nail Nipper
Curved Jaw
Double Spring
Size: 11.5cm



LSI-00-236
Ingrown Diabetic
Nail Nipper
Double Spring
Size: 11.5cm



LSI-00-237
Ingrown Nail Nipper
Extra Narrow Jaw
Double Spring
Size: 11.5cm



LSI-00-238
Ingrown Nail Nipper
Concave Jaw
Double Spring
Size: 11.5cm




Ingrown Nail Nippers



Box-joint





LSI-00-239
Ingrown Nail Nipper
Standard Jaw
Double Spring
Size: 13cm




LSI-00-240
Ingrown Nail Nipper
Pointed Jaw
Double Spring
Size: 13cm



LSI-00-241
Ingrown Nail Nipper
Slim Jaw
Double Spring
Size: 13cm




LSI-00-242
Ingrown Nail Splitter
Double Spring
Size: 13cm






LSI-00-243
Ingrown Nail Nipper
Curved Jaw
Double Spring
Size: 13cm



LSI-00-244
Ingrown Diabetic
Nail Nipper
Double Spring
Size: 13cm



INGROWN NAIL NIPPER

EXTRA NARROW JAW



LSI-00-245
Ingrown Nail Nipper
Extra Narrow Jaw
Double Spring
Size: 13cm



Box-joint




Ingrown Nail Nippers




Box-joint



LSI-00-246
Ingrown Nail Nipper
Standard Jaw
Double Spring
Size: 15cm



LSI-00-247
Ingrown Nail Nipper
Slim Jaw
Double Spring
Size: 15cm



LSI-00-248
Ingrown Nail Nipper
Pointed Jaw
Double Spring
Size: 15cm




LSI-00-249
Ingrown Nail
Curved Jaw
Double Spring
Size: 15cm






PODOLOGY SET





LSI-00-260
Ingrown Nipper
Scissor Ring
Size: 10cm



LSI-00-262
Ingrown Nipper
Scissor Ring
Size: 8.5cm



LSI-00-263
Flat Jaw Nipper
Scissor Ring
Size: 10cm



LSI-00-264
Flat Jaw Nipper
Scissor Ring
Size: 8.5cm





LSI-00-265
Nail Clipper
Curved
Multicolor Finish
Size: 7cm




LSI-00-266
Nail Clipper
Curved
Satin Finish
Size: 7cm




LSI-00-267
Nail Clipper
Straight
Satin Finish
Size: 7cm

Comfortable and Effortless
Feel The Quality




Nail Clippers





LSI-00-268
Conical Flat Nose Plier With Serration

wire up to .028"
(0.7mm.) max.




LSI-00-269
Long Beak Light Wire Plier

wire up to .020"
(0.5mm.) max.



LSI-00-270
Universal General Use Plier





Both Tips are serrated. Bends hardwire
up to .039" (1.0mm) max.
Cuts hardwire up to .028" (0.7mm)
max.





LSI-00-271
Jarabak Plier with Cutter

wire up to .020"
(0.5mm.) max..



WE BRING THEM TO LIFE



Box-joint






LSI-00-272

Hard Wire Cutter Straight Jaws

wire up to .028"
(0.7mm.) max.



**Heavy Duty
Power Cutter**



LSI-00-273

Hard Wire Cutter-Lab Tool

Max .050"



German



CUTICLE NIPPERS


LASANI SURGICAL Cuticle nippers are manufactured using high-pressure forging of layers made of **100%** surgical grade stainless steel. Our nippers are designed for professionals to use for durability, sharp cutting and long-lasting sharpness. The quality of the steel has a crucial influence on the quality of the nippers.

Our Cuticle nippers are box-jointed and have the quality to stop side-way movements and adjusts the alignment of the cutting edges. The handles have a round tension to give you the highest possible grip which allows the hand to have maximum control.


Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, it features gently curved, palm-friendly, easy-to-grip ergonomically designed satin finish which simplifies easy cleaning for optimal comfort, performance and expresses the ultimate in functional beauty.



Total Length 11.5cm



Total Length 11.5cm



Total Length 10cm



Double Spring



Single Spring



Lasani Surgical Cuticle nippers come in **4 unique handle** sizes and make the nippers user friendly.



Box-joint

All of Lasani Surgical nippers have springs that are ass embled with a screw which can be replaced.



LSI-00-274

Cuticle Nipper 6mm Jaw

Double Spring
Size: 10cm



LSI-00-275

Cuticle Nipper 4mm Jaw

Double Spring
Size: 10cm



LSI-00-276

Cuticle Nipper 3mm Jaw

Double Spring
Size: 10cm



LSI-00-277

Cuticle Nipper 7mm Jaw

Double Spring
Size: 10cm



LSI-00-278

Cuticle Nipper 9mm Jaw

Double Spring
Size: 10cm





LSI-00-283

Cuticle Nipper 6mm Jaw

Single Spring

Size: 10cm



LSI-00-279

Cuticle Nipper 4mm Jaw

Single Spring

Size: 10cm



LSI-00-280

Cuticle Nipper 3mm Jaw

Single Spring

Size: 10cm



LSI-00-281

Cuticle Nipper 7mm Jaw

Single Spring

Size: 10cm



LSI-00-282

Cuticle Nipper 9mm Jaw


Single Spring

Size: 10cm

Cuticle Nippers



Box-joint



LSI-00-284

Cuticle Nipper 6mm Jaw

Double Spring
Size: 11.5cm



LSI-00-285

Cuticle Nipper 4mm Jaw

Double Spring
Size: 11.5cm



LSI-00-286

Cuticle Nipper 3mm Jaw

Double Spring
Size: 11.5cm



LSI-00-287

Cuticle Nipper 7mm Jaw

Double Spring
Size: 11.5cm



LSI-00-288

Cuticle Nipper 9mm Jaw

Double Spring
Size: 11.5cm





LSI-00-289

Cuticle Nipper 6mm Jaw

Single Spring

Size: 11.5cm



LSI-00-290

Cuticle Nipper 4mm Jaw

Single Spring

Size: 11.5cm



LSI-00-291

Cuticle Nipper 3mm Jaw

Single Spring

Size: 11.5cm



LSI-00-292

Cuticle Nipper 7mm Jaw

Single Spring

Size: 11.5cm



LSI-00-293

Cuticle Nipper 9mm Jaw

Single Spring

Size: 11.5cm





LSI-00-294
Cuticle Nipper Fino
8mm Jaw
Double Spring
Size: 11.5cm



LSI-00-295
Cuticle Nipper Fino
6mm Jaw
Double Spring
Size: 11.5cm



LSI-00-296
Cuticle Nipper Fino
4mm Jaw
Double Spring
Size: 11.5cm



LSI-00-297
Cuticle Nipper Fino
3mm Jaw
Double Spring
Size: 11.5cm



LSI-00-298
Cuticle Nipper Fino
8mm Jaw
Single Spring
Size: 11.5cm



LSI-00-299
Cuticle Nipper Fino
6mm Jaw
Single Spring
Size: 11.5cm



LSI-00-300
Cuticle Nipper Fino
4mm Jaw
Single Spring
Size: 11.5cm



LSI-00-301
Cuticle Nipper Fino
3mm Jaw
Single Spring
Size: 11.5cm



Cuticle Nippers



Box-joint



LSI-00-302
Cuticle Nipper
8mm Jaw
Double Spring
Size: 11.5cm



LSI-00-303
Cuticle Nipper
6mm Jaw
Double Spring
Size: 11.5cm



LSI-00-304
Cuticle Nipper
3mm Jaw
Double Spring
Size: 11.5cm



» INNOVATION
» CLASS
» QUALITY





LSI-00-305

Acrylic Nail Nipper

Single Spring
Jaw Size: 6mm
Size: 10.5cm




Box-joint





LSI-00-306
Tip Slicer Straight Blade
Size: 13cm



Acrylic Nail Set



REF: L-6170



LSI-00-314

LSI-00-313

LSI-00-312

LSI-00-311


LSI-00-310



LSI-00-307
10.5cm



LSI-00-308
11.5cm



LSI-00-309
11.5cm with 7mm Jaw



LSI-00-315
Cuticle Nipper
Scissor Ring

Size: 10cm
jaw size: 5mm

Also Available in 7mm jaw size
REF: L-6173 jaw size: 7mm



LSI-00-316
Cuticle Nipper
Scissor Ring

Size: 8.5cm
jaw size: 6mm



LSI-00-317
Cuticle Nipper
12mm Jaw
Double Spring
Size: 12cm




LSI-00-318
Cuticle Nipper
15mm Jaw
Double Spring
Size: 13cm




Cuticle Nippers



Box-joint



LSI-00-319
Cuticle Nipper
6mm Jaw
Single Spring
Size: 10cm



LSI-00-320
Cuticle Nipper
4mm Jaw
Single Spring
Size: 10cm




LSI-00-321
Cuticle Nipper
8mm Jaw
Single Spring
Size: 10cm



Slant blade
will not block the view

*blade length 6mm



CUTICLE PUSHERS

LASANI SURGICAL Cuticle Pushers are designed to safely remove dead skin from the nail bed without causing damage or scratches to the nail. They allow cuticles to be gently pushed back, helping nails grow stronger and healthier.

LASANI SURGICAL Cuticle Pushers come in a wide variety of forms and models. These include small to medium scoops, ingrown nail lifters, pushers with a knife edge, flat-edge pushers, extractors, arrow-point tips, and pterygium removers. LASANI SURGICAL tools are manufactured in diverse model ranges to ensure comfort and user-friendliness for their intended purposes.

Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, our pushers feature many types of edges, palm-friendly, easy-to-grip ergonomically designed finishes on the handles which gives optimal comfort, performance and expresses the ultimate in functional beauty. Our non-slip grooved texture on the handles allow a comfortable grip while the tool is in use.





Long Jaw

LSI-00-322

Cuticle Pusher with Ingrown Lifter Long Jaw



LSI-00-323

Cuticle Pusher with Angled Sharp Edge



LSI-00-324

Cuticle Pusher with Flat Edge



LSI-00-325

Cuticle Pusher Narrow with Flat Edge



LSI-00-326

Cuticle Pusher with Small Curved Pterygium Remover



LSI-00-327

Cuticle Pusher Rounded and Bent Blade



LSI-00-328

Cuticle Pusher with Nail Edge Cleaner 2mm



LSI-00-329
Cuticle Flat Pusher and Bent
Blade



LSI-00-330
Cuticle Pusher with Knife



LSI-00-331
Cuticle Pusher Small Scoop
with Knife



LSI-00-332
Cuticle Pusher Small Scoop
with Knife



LSI-00-333
Cuticle Pusher Small Scoop and Remover



LSI-00-334
Cuticle Pusher with Flat End



LSI-00-335
Cuticle Pusher Small Scoop with Knife

Cuticle Pushers



LSI-00-336
Flat Edges Cuticle Pushers




LSI-00-337
Cuticle Pusher Double Ended



LSI-00-338
Arrow Pointed Tip with
Flat Pusher



LSI-00-339
Cuticle Pusher Angled Scoop
with Knife



LSI-00-340
Ingrown Nail Lifter and Angled Scoop



LSI-00-341
Triangle Cuticle Peeler Scraper



LSI-00-342
Cuticle Pusher with Small Curved Pterygium Remover Left

EXCAVATORS, CURETTES AND PROBES



LASANI SURGICAL Excavators are designed to safely remove dead skin from the side walls of ingrown nails.

LASANI SURGICAL Curettes are designed to clean under fingernails and toenails and are also used to remove cuticle from the nail plate.

LASANI SURGICAL Excavators and Curettes are available in a wide variety of forms and model ranges. Excavators come in single- or double-sided designs, with round edges and various sizes. Our curettes include Volkman curettes, curettes without holes, oval-shaped curettes, and more. LASANI SURGICAL tools are manufactured in diverse model ranges to ensure comfort and user-friendliness for their intended purposes.

Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, our excavators and curettes feature many types of edges, palm-friendly, easy-to-grip ergonomically designed finishes on the handles which gives optimal comfort, performance and expresses the ultimate in functional beauty. Our non-slip grooved texture on the handles allow a comfortable grip while the tool is in use.

Available in Different Handles





LSI-00-343
Nail Cleaner Sharp Edges Curved



LSI-00-344
Nail Curette Angled 2 mm
& 3 mm, 13cm



LSI-00-345
Nail Excavator Double Sided



LSI-00-346
Nail Excavator Round Edge
D/Ended 13cm



LSI-00-347
Nail Pusher with Knife double ended



LSI-00-348
Chisel with Knife Double Ended



LSI-00-349
Nail Curette with Flat Square Edge 13cm



LSI-00-350
Curette Excavator, Double End,
with Holes, 3 mm and 4 mm




LSI-00-351
Curette Excavator, Double End,
with Holes, 1.5 mm and 2.5 mm



LSI-00-352
Curette Excavator, Double End,
with Holes, 1 mm and 2 mm




LSI-00-353
Nail Curette Angled 2 mm
& 3 mm, 13cm



LSI-00-354
Nail Probe Angled Flat fine Head 14cm



LSI-00-355
Nail Probe Swan Neck Flat 14cm



LSI-00-356
Nail Probe Sharp 14cm



LSI-00-357
Nail Probe & Excavator

Excavators & Probes



LSI-00-358
Chisel Sharp Single sided
Slightly Curved



LSI-00-359
Chisel Sharp Single sided
Angled



LSI-00-360
Chisel Sharp Fine Single
Sided 2 mm



LSI-00-361
Chisel Sharp Fine Single
Sided 1 mm



LSI-00-362
Nail Excavator Round Edge D/Ended 13cm



LSI-00-363
Nail Excavator and Arrow Point Sharp



LSI-00-364
Nail Excavator Double Sided

EXTRACTORS

LASANI SURGICAL Extractors are designed to safely remove blackheads and whiteheads from the skin. LASANI SURGICAL Extractors are available in a wide variety of forms and model ranges. These include double-sided blunt hole extractors, needle extractors, blackhead extractors, round scoop extractors, lancet extractors, comedone extractors, and more. LASANI SURGICAL tools are manufactured in diverse model ranges to ensure comfort and user-friendliness for their intended purposes.

Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, our extractors feature many types of edges, palm-friendly, easy-to-grip ergonomically designed finishes on the handles which gives optimal comfort, performance and expresses the ultimate in functional beauty. Our non-slip grooved texture on the handles allow a comfortable grip while the tool is in use.





LSI-00-365
Extractor Lancet Double
Sided Removable Head



LSI-00-366
Extractor with 90° Sharp Needle



LSI-00-367
Extractor Double Sided



LSI-00-368
Angled Extractor with
Langenbeck Curette Large



LSI-00-369
Angled Extractor Double Sided 1.5mm & 2.5mm



LSI-00-370
Extractor Double Sided 1.5mm & 2.5mm



LSI-00-371
Langenbeck Curette Large



LSI-00-372
Round Scoop Double
Ended Extractor



LSI-00-373
Oval Scoop Extractor with
Langenbeck Curette Large



LSI-00-374
Angled with Large Oval
Scoop Extractor



LSI-00-375
Extractor Double Sided For
Bulk & Single Black Heads



LSI-00-376
Straight & Angled Extractor Double Sided For Bulk & Single Black Head



LSI-00-377
Comedone Extractor Curved Blackhead Tweezer Stainless Steel Size; 6"



LSI-00-378
Chisel Handle Round for
Microsurgery Blade Mini




LSI-00-379
Chisel Handle Round for
Microsurgery Blade



LSI-00-380
Chisel Handle Round for
Microsurgery Blade No. 2



LSI-00-381
Scalpel Handle No. 4



LSI-00-382
Scalpel Handle No. 3






LSI-00-383
Scalpel Handle No. 7

INGROWN NAIL FILES

Nothing feels as satisfying as a good nail file, and LASANI SURGICAL high-tech nail files far surpass expectations. Manufactured in many unique shapes, usage methods, and surfaces, these exquisitely crafted files feature one abrasive side for precise filing and one smoother side for finishing touches. When held in hand, you'll notice a finger dip that creates a comfortable, ergonomic grip.

LASANI SURGICAL Nail Files come in a wide variety of forms and model ranges. These include diamond dust files, finger and toenail files, glass nail files, ingrown nail files, straight-head fine-cut files, black nail files, angled nail files, ultra-fine nail files, and many more. All are designed to maximize comfort and ensure user-friendliness for their intended purposes.

Available Different Handles





LSI-00-384
Nail Rasp Backward Cutting
Small Head and Curved Head



LSI-00-385
Black Ingrown Nail Large
Head straight and Angle



LSI-00-386
8 Black Ingrown Nail File Oval
Head Straight and Angled



LSI-00-387
Nail File Small Fine Head &
Excavator Hollow Handle 15cm



LSI-00-388
Nail File backward Cutting Curved Hollow Handle 16cm



LSI-00-389
Nail File broad Fine Cut STR & CVD Head Hollow Handle 15cm



LSI-00-390
Nail File Supper Fine Head Straight and Angled Head Hollow Handle 15cm



LSI-00-391
Nail File Large Head CVD
Course Cut 14cm



LSI-00-392
Nail File Small Oval
Head Fine Cut



LSI-00-393
Nail File Medium Oval
Head Fine Cut



LSI-00-394
Nail File Large Oval
Head Fine Cut



LSI-00-395
Nail File Supper Fine Head Fine Cut STR 14cm



Ingrown Nail Files



LSI-00-396
Arrow Head Nail File Curved
& Straight



LSI-00-397
Oval Head Nail File Curved
& Straight



LSI-00-398
Oval Head Nail File with
Curette Round



LSI-00-399
Ingrown Narrow Head Nail File
Straight Coarse & Fine



LSI-00-400
Ingrown Narrow Head Nail File Fine Curved & Straight



LSI-00-401
Standard Head Nail File Coarse Curved & Straight



LSI-00-402
Angled Head Nail File Fine



LSI-00-403
Straight Head Nail File with
Sharp Curette



LSI-00-404
File with Nail Spatula



LSI-00-405
Ingrown Straight Nail File Double
Sided Fine with Curette Round



LSI-00-406
Angled Head Nail File Coarse



LSI-00-407
Straight Serrated Nail File Angled



LSI-00-408
Cross Serrated Nail File Angled

Ingrown Nail Files



LSI-00-409
Ingrown Nail File Fine Serrated
Angled Head




LSI-00-410
Ingrown Nail File Coarse
Single Sided Straight Small




LSI-00-411
Ingrown Nail File Coarse
Single Sided Straight Medium




LSI-00-412
Ingrown Nail File Coarse
Single Sided Straight Large



LSI-00-413
Ingrown Nail File Fine Serrated Narrow Head Small



LSI-00-414
Ingrown Nail File Fine Serrated Narrow Head Medium



LSI-00-415
Ingrown Nail File Fine Serrated Narrow Head Large



LSI-00-416
Ingrown Nail File Fine Serrated
Arrow Head Curved




LSI-00-417
Ingrown Nail File Fine Serrated
Arrow Head Straight Large




LSI-00-418
Ingrown Nail File Fine Serrated
Arrow Head Straight Medium



LSI-00-419
Ingrown Nail File Fine Serrated
Arrow Head Straight Small



LSI-00-420
Oval Head Nail File Straight Small Single Sided



LSI-00-421
Oval Head Nail File Straight Medium Single Sided



LSI-00-422
Oval Head Nail File Straight Large Single Sided

Ingrown Nail Files



SPATULAS & NAIL PACKERS

Lasani Surgical Spatulas are designed to spread a thin layer of warmed wax onto your skin, in the direction of hair growth.

LASANI SURGICAL Spatulas are available in a variety of forms, including round circle spatulas, slim spatulas, wide-blade spatulas, double-sided waxing spatulas, flat curved spatulas, wooden-handle spatulas, and stainless-steel handle spatulas. LASANI SURGICAL tools are manufactured in diverse model ranges to ensure comfort and user-friendliness for their intended purposes.

Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, our spatulas feature many types of edges, palm-friendly, easy-to-grip ergonomically designed finishes on the handles which gives optimal comfort, performance and expresses the ultimate in functional beauty. Our non-slip grooved texture on the handles allow a comfortable grip while the tool is in use.

Available in 2 Different Handles





LSI-00-423
Nail Spatula Hollow Handle
6mm Broad Blade



LSI-00-424
Nail Spatula Hollow Handle
3mm Broad Blade



LSI-00-425
Elevator & Nail Packer



LSI-00-426
Nail Spatula Double Sided



LSI-00-427
Elevator & Nail Packer



LSI-00-428
Waxing Spatula Double Sided



LSI-00-429
Nail Spatula with Flat Probe 15cm



LSI-00-430
Nail Spatula with Slim Flat Probe

Spatulas & Nail Packers



LSI-00-431
Excavator with 2X25mm
Packer Blunt Hook Spatula



LSI-00-432
Spatula & Nail Packer




LSI-00-433
Arrow Point & Round Circle
Shape Spatula



LSI-00-434
Straight Probe and Blunt Hook
Size; 14cm



LSI-00-435
Nail Spatula 2mm Single Ended 13cm



LSI-00-436
Nail Spatula 4mm Single Ended 13cm




LSI-00-437
Nail Spatula 6mm Single Ended 13cm




LSI-00-438
Spatula Wood Handle



LSI-00-439
Stainless Steel Spatula
Comfortable Grip



LSI-00-440
Waxing Spatula Single Sided



LSI-00-441
Mirror 24mm #5
with Solid Handle

Spatulas



EYEBROW & EYELASH TWEEZERS

Lasani Surgical Tweezers are manufactured using high-pressure forging of layers made of 100% surgical-grade stainless steel, which is capable of performing under normal or expected conditions. Our tweezers are designed for professionals to ensure durability and long-lasting ease of tension. The quality of the steel has a crucial influence on the overall quality of the tweezers.

Lasani Surgical Tweezers come in various forms with a diverse model range. This includes slanted tweezers, narrow slanted tweezers, straight claw tweezers, wide-handle tweezers, electrolysis tweezers in various sizes, body waxing tweezers, eyelash pointed angled tweezers, ultra-sharp pointed curved tweezers, angled needle-nose tweezers, eyelash extension tweezers, and many more. Lasani Surgical tools are manufactured in diverse model ranges to provide comfort and user-friendly performance for their intended purpose.

Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, our tweezers feature many types of edges, palm-friendly, easy-to-grip ergonomically designed finishes on the handles which gives optimal comfort, performance and expresses the ultimate in functional beauty. Our tweezer joints are perfectly jointed on the handles to allow a comfortable grip while the tool is in use and allow a decent amount of tension while in use.



Standard Tip



LSI-00-442
Eyebrow Tweezer
Slanted Tip
Size: 10cm



LSI-00-443
Eyebrow Tweezer
Slanted Tip with Hole Grip
Size: 10cm



LSI-00-444
Eyebrow Tweezer
Straight
Size: 10cm



LSI-00-445
Eyebrow Tweezer
Slanted with Wide Grip
Size: 10cm



LSI-00-446
Eyebrow Tweezer
Angled with Comb
Size: 10.5cm



LSI-00-447
Eyebrow Tweezer
Curved
Size: 10cm

Eyebrow Tweezers Standard



LSI-00-448
Eyebrow Tweezer
Slanted Angled
Size: 10cm




LSI-00-449
Eyebrow Tweezer
Slanted with Oval Grip
Size: 10cm



LSI-00-450
Eyebrow Tweezer Straight
with Standard Grip
Size: 10.5cm

Standard Tip

EYEBROW TWEEZER DRAWING



Standard Tip



LSI-00-451
Eyebrow Tweezer
Slanted
Size: 10cm



LSI-00-452
Eyebrow Tweezer
Straight Narrow
Size: 9cm



LSI-00-453
Eyebrow Tweezer
Slanted
Size: 10cm



LSI-00-454
Eyebrow Tweezer
Slanted with Hole Grip
Size: 10cm



LSI-00-455
Eyebrow Tweezer
Slanted with Wide Grip
Size: 10cm



Eyebrow Tweezers



Narrow Tip



LSI-00-456
Eyebrow Tweezer
Slanted Narrow
Size: 9cm



LSI-00-457
Eyebrow Tweezer Slanted
Narrow with Standard Grip
Size: 9cm



LSI-00-458
Eyebrow Tweezer
Straight Narrow
Size: 9cm



LSI-00-459
Eyebrow Tweezer Slanted
Narrow with Cone Grip
Size: 10cm



LSI-00-460
Eyebrow Tweezer
POINTED
Size: 9cm



LSI-00-461
Body Waxing Tweezer
Size: 9cm





LSI-00-462
Eyelash Tweezer Straight
Narrow Pointed
Size: 11.5cm



LSI-00-463
Eyelash Tweezer
Straight Fine Pointed
Size: 11.5cm



LSI-00-464
Eyelash Tweezer
Angled Fine Pointed
Size: 11.5cm



LSI-00-465
Eyelash Tweezer Double
Curved Fine Pointed
Size: 13cm



LSI-00-466
Eyelash Tweezer Angled
45° Degree
Size: 11.5cm



LSI-00-467
Eyelash Tweezer Angled
90° Degree
Size: 11.5cm



LSI-00-468
Eyelash Tweezer Extra
Fine Slightly Angled
 Size: 11.5cm



LSI-00-469
Eyelash Tweezer Extra
Fine Curved
 Size: 11.5cm



LSI-00-470
Eyelash Tweezer Extra
Fine Curved
 Size: 14cm





LSI-00-471
Eyelash Tweezer Extra
Fine Angled
 Size: 11.5cm



LSI-00-472
Eyelash Tweezer Extra
Fine Slightly Angled
 Size: 11.5cm



LSI-00-473
Eyelash Tweezer
Straight Fine Pointed
 Size: 11.5cm






LSI-00-479
Adson Tissue Forceps
 Available in two Sizes
 12cm & 15cm



LSI-00-480
Adson Tissue Forceps
 Size: 12cm



LSI-00-481
Desecting Fine Tissue Forceps Curved
 Size: 11.5cm



LSI-00-482
Dissecting Waugh Tissue Forceps
 Size: 13cm



LSI-00-483
Dressing Tissue Forceps Straight
 Size: 14cm



LSI-00-484
Adson Tissue Forceps Arrow Tips
 Size: 12cm



NAIL. CUTICLE AND EYEBROW SCISSORS

Lasani Surgical Scissors are designed to safely remove hair from the body and/or cuticles from the nail bed without causing damage or scratches to the nail. They allow precise cuticle trimming, helping promote stronger and healthier nail growth.

Lasani Surgical Scissors come in various forms with a diverse range of models. This includes stork scissors, curved nail scissors, toenail scissors, fiberglass scissors, ultra-fine scissors, mustache scissors, safety scissors, fine and curved pointed cuticle scissors, and more. Lasani Surgical tools are manufactured in diverse model ranges to ensure comfort and user-friendly performance for their intended purpose.

Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, our scissors feature many types of blades, palm-friendly, easy-to-grip ergonomically designed finishes on the rings which gives optimal comfort, performance and expresses the ultimate in functional beauty. Our sand finished texture on the handles allow a comfortable grip while the tool is in use.





LSI-00-485
Arrow Pointed Scissor
Curved
Size: 9cm



LSI-00-486
Arrow Pointed Scissor
Double Curved
Size: 9cm



LSI-00-487
Narrow Point Nail
Scissor Curved
Size: 9cm



LSI-00-488
Fine Scissor Curved
Size: 9cm



LSI-00-489
Nail Scissor Curved
Size: 9cm



LSI-00-490
Fine Scissor Curved
Size: 10cm



LSI-00-491
Embroidery and Sewing
Scissor
Size: 9cm



LSI-00-492
Nose Hair Scissor
Size: 11.5cm



LSI-00-493
Dissecting Scissor
Size: 11.5cm



LSI-00-494
Large Rings Embroidery
Scissor
Size: 10cm



LSI-00-495
Embroidery Scissor
Size: 9cm



LSI-00-496
Embroidery Scissor
Straight
Size: 10cm



LSI-00-497
Stork Scissor Arrow
Pointed
Size: 9cm



LSI-00-498
Stork Scissor Arrow
Pointed Long Blade
Size: 11.5cm




LSI-00-499
Embroidery Scissor
Straight
Size: 9cm



LSI-00-500
Micro Spring Scissor
Straight
Size: 12cm



LSI-00-501
Micro Spring Scissor
Curved
Size: 12cm



Sharp Blades



LSI-00-502
Ultra Fine Scissors
Straight
Size: 9cm



LSI-00-503
Ultra Fine Scissors
Straight
Size: 10cm



LSI-00-504
Ultra Fine Scissors
Straight
Size: 11.5cm



LSI-00-505
Ultra Fine Scissors
Curved
Size: 9cm



LSI-00-506
Ultra Fine Scissors
Curved
Size: 10cm



LSI-00-507
Ultra Fine Scissors
Curved
Size: 11.5cm



LSI-00-508

**Super Operating Scissor
Straight Sharp/Sharp**

Size: 15cm



LSI-00-509

**Operating Scissor Straight
Sharp/Blunt**

Size: 15cm



LSI-00-510

**Operating Scissor Straight
Blunt/Blunt**

Size: 15cm





LSI-00-511
Supper Edge Iris Scissor
Curved
Size: 11.5cm



LSI-00-512
Spencer Stitch Scissor
Curved Delicate
Size: 10.5cm



LSI-00-513
Arrow Pointed Scissor
Double Curved
Size: 11.5cm





LSI-00-514
Narrow Point Nail
Scissor Double Curved
Size: 12cm




LSI-00-515
Lister Bandage Scissor
Size: 15cm










LSI-00-519
Mosquito Forceps
Size: 12.5cm



LSI-00-520
Mosquito Forceps
Size: 10cm



Standard
Tip

LSI-00-521
Mathieu Standard Tip
Double Spring


Mosquito Forceps



STAINLESS STEEL FOOT FILE

Looking to invest in a foot file for long-term use without worrying about corrosion? Lasani Surgical currently offers two highly demanded stainless-steel foot files. Our Lasani Surgical foot file is designed to accommodate coarse and medium grit adhesive stickers, making it ideal for use on multiple clients within a short span of time. The file can be purchased with coarse or medium grit stick-and-peel sandpapers, available in packs of 50 each, allowing up to 50 customers to be serviced in a spa or academy using just one tool. Another highly popular option is the Lasani Surgical solid stainless-steel foot file, manufactured using high-pressure forging of layered 100% surgical-grade stainless steel. This file is capable of performing under normal or expected professional conditions and is designed for durability and long-lasting ease of use. The quality of the steel plays a crucial role in ensuring the superior performance and longevity of the foot file.

Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, our stainless steel file features palm-friendly, easy-to-grip ergonomically designed finishes on the handles which gives optimal comfort, performance and expresses the ultimate in functional beauty. This file can be auto claved, reused and is durable for longevity for thousands of uses without any fear of wear or damage,



LSI-00-522
Stainless Steel Foot File
Straight





STAINLESS STEEL FOOT FILE

Looking to invest in a foot file for long-term use without worrying about corrosion? Lasani Surgical currently offers two highly demanded stainless-steel foot files. Our Lasani Surgical foot file is designed to accommodate coarse and medium grit adhesive stickers, making it ideal for use on multiple clients within a short span of time. The file can be purchased with coarse or medium grit stick-and-peel sandpapers, available in packs of 50 each, allowing up to 50 customers to be serviced in a spa or academy using just one tool. Another highly popular option is the Lasani Surgical solid stainless-steel foot file, manufactured using high-pressure forging of layered 100% surgical-grade stainless steel. This file is capable of performing under normal or expected professional conditions and is designed for durability and long-lasting ease of use. The quality of the steel plays a crucial role in ensuring the superior performance and longevity of the foot file.

Our highly redefined passivation process protects instruments against corrosion due to repeated sterilization. Responding to customer's needs, our stainless steel file features palm-friendly, easy-to-grip ergonomically designed finishes on the handles which gives optimal comfort, performance and expresses the ultimate in functional beauty. This file can be auto claved, reused and is durable for longevity for thousands of uses without any fear of wear or damage,



LSI-00-523
Stainless Steel Foot File
Straight



Coarse



Medium



WOODEN FOOT FILE


Lasani Surgical is very proud of its top-selling environmentally friendly wooden foot file. It is designed to provide a comfortable grip during use, featuring tough and durable filing surfaces for effective foot care. The file consists of two sides: one side for removing calluses and the other for smoothing and buffing the skin. The construction of our wooden foot file undergoes approximately 50 manufacturing processes, making it waterproof and easily washable without fear of damage. Our wooden files are fully sanitizable, helping resist bacteria and other microbes on the soles of the feet. Advanced manufacturing technology ensures durability, with multiple processes reinforcing the longevity of the file. The rounded edges provide added comfort for users, whether for personal home use or professional spa treatments. When used dry, the file effectively removes calluses; when used wet, it helps buff and smooth the skin.

LSI-00-524


Wooden Foot File Straight

Waterproof Gauranteed
Tender Heel
Remove Dead Skin Softly And Smoothly
Remove Calluses
Fine And Smooth, No Hurt To The Feet





LSI-00-525
Wooden Foot File
Thick Handle
Size: 11.5"



LSI-00-526
Plastic Handle Foot
File Straight
Size: 14cm





REF: L-6584
Callus Shaver
Stainless Steel



insert now blade.
Hold sides and press
down on handle and
Pull towards You.



Hold Sides Push
Handle forward



Lay on Flat
Surface

LSI-00-527
**Foot Rasps 4 Different
Surface**
8.5mm Broad




LSI-00-528
Foot Rasps Single Sided
Fine & Coarse



LSI-00-529
Foot Rasps Single Sided
Fine



LSI-00-530
Foot Rasps Peel & Stick
Stainless Steel File



LSI-00-531
Foot File Disposable Paper
100 Grit (Medium)
Pack of 50 Pieces




LSI-00-532
Foot File Disposable Paper
100 Grit (Coarse)
Pack of 50 Pieces



LSI-00-533
Foot Rasps Peel & Stick
Stainless Steel File



Coarse



Medium

LSI-00-534
Rifflor Angled Foot
Rasps Size: 20cm



LSI-00-535
Diamond Deb Foot
Dresser File
 Size: 20cm



LSI-00-536
Diamond Deb Foot
Dresser File
 Size: 15cm



LSI-00-537
Serrated Foot Dresser
File
 Size: 15cm



LSI-00-538
Nail File Double Sided
Stainless Steel Coarse
/Medium



LSI-00-539
Nail File Double Sided
Stainless Steel Coarse/Medium



LSI-00-540
Nail File Single Sided
Acrylic



LSI-00-541
Diamond Deb Nail File
Double Sided Stainless
Steel Coarse/Medium
Size: 18cm




LSI-00-542
Diamond Deb Nail File
Double Sided Stainless
Steel Coarse/Medium
Size: 15cm



LSI-00-543
Nail File Double Sided
Stainless Steel Coarse
/Medium
Size: 15cm




LSI-00-544
Diamond Deb Nail File
Double Sided Stainless
Steel Size: 10cm



LSI-00-545
Serrated Nail File
 Size: 10cm



LSI-00-546
Serrated Nail File
 Size: 15cm



LSI-00-547
Serrated Nail File
 Size: 10cm



LSI-00-548
**Diamond Deb Foot
 Dresser Slim**
 Size: 20cm



LSI-00-549
**Diamond Deb Foot
 Dresser Slim 10cm**
 Size: 15cm



LSI-00-550
**Diamond Deb Foot
 Dresser Slim 10cm**
 Size: 10cm





LSI-00-551

Kidney Tray Large

Size: 30.48 x 12.19 x 5.08cm



LSI-00-552

Kidney Tray Medium

Size: 26.67 x 15.24 x 4.45cm



LSI-00-553

Kidney Tray Small

Size: 21.59 x 12.07 x 3.81cm



LSI-00-554

Kidney Tray Extra Small

Size: 15.88 x 8.89 x 2.54cm



LSI-00-555

Kidney Tray Small

Size: 20.96 x 11.81 x 1.65cm



Available in two Sizes


LSI-00-556

Instruments Tray

Size: 30.48 x 20.32 x 1.27cm

Instruments Tray

Size: 33x23x5cm



LSI-00-557
Sterilization Trays
Size: 25X25X5 Inches

LSI-00-558
Sterilization Trays
Size: 10X5X2.25 Inches



LSI-00-559
Sterilization Trays
Size: 8X4X2 Inches



LSI-00-560
Sterilization Trays
Size: 12X6X2.50 Inches



LSI-00-561
**Podology Set Black/
Black**

Inside Set:
11 Instruments





LSI-00-562

**Podology Set Camel/
Camel**

Inside Set:

12 Instruments



LSI-00-563

**Podology Set Black/
Black**

Inside Set:

6 Instruments





LSI-00-564
Pocket manicure set
Camel/Black
Inside Set:
6 Instruments





SI-00-565

Pocket manicure set
Black/Black

Inside Set:
Instruments



LSI-00-566

Pocket manicure set
Black/Black


Inside Set:
4 Instruments

Notes




Lasani Surgical

 Sialkot, Punjab, Pakistan

 (052) 355-87-36

 Info@lasanisurgicalind.com

 www.lasanisurgicalind.com

