

Since 1998 we have served the HealthCare, Surgical Industry with the finest products and saved the lives of people. Lasani Surgical has proven itself as a quality manufacturer and exporter of Dental, Surgical & Veterinary instruments

Lasani Surgical aims to contribute and help the world of health care, with the most innovative technology in the Health Care Industry. Lasani Surgical is aiming to be a solid global company with high class quality with a passion for manufacturing and spirit for meeting a challenge and always looking for new innovations. Our team is always in contact with the leading doctors from the globe to work on universal ideas.

We stand behind every product we offer, listen to our customers and exceed their expectations. We are a rapidly growing company, achieving high Quality standards in a very short period. Lasani Surgical is a truly customer centric company and focus on creation of new values.



Introduction

Quality Policy:

-Everybody talks about Quality but Lasani Surgical provides Quality. We consider the details. Locks, catches and box lock mechanisms are electro-polished and Forceps and Tweezers are laser-welded at the back.

-We proudly mark our instruments with a limited lifetime warranty from rust and craftsmanship.

-All the instruments are made using AISI 420, 410 & 18/18 (304) grade stainless steel. We strictly do not manufacture our instruments with economy grade Steel. We have been importing the raw material from the world's renown companies since 1975. All the instruments are manufactured with Japanese and German stainless steel and finalised with German finishing.

-The Directors, management & staff makes sure that all the products are consistently checked and reviewed at all the stages to achieve the quality mandatory for our customers' satisfaction.

Lasani Surgical Production Standards:

Lasani Surgical applied international standards during the process of manufacturing. Following a list and reference of those standards for our customers.

- Acquiring quality stainless steel as per ASTM standards for our customers.
- Acquiring quality stainless steel as per ASTM standard F899-95 & ISO 713501.
- Corrosion of Instruments in conformity of ASTM Standard F1089-02
- Passivation in Conformity to ASTM Standard B967-05.
- Electro polishing as per ASTM Standard B912-02
- Rockwell Hardness test in accordance with ASTM standard E18-00.
- Cleaning and Descaling of stainless steel as per ASTM standard A380-99 (R05)




THE SIALKOT CHAMBER OF COMMERCE & INDUSTRY	
Member Detail	
Contact Person:	NOUMAN ALI
Company Name:	LASANI SURGICAL INDUSTRIES
Member Class:	Associate
Member ID:	A-33124
Designation:	Partnership
Category:	GENERAL GOODS
Address:	JOURIAN KALAN, P.O. KHAS, WAZIRABAD ROAD, SIALKOT
Mobile:	0300-6141183
Email:	lasurgical@gmail.com, lasurgical@msl.com
Word:	

LOCKING PLATES	2-16
LOCKING SCREW	17-18
BONE PLATE	19-30
BONE SCREW	31-39
INTERLOCKING NAIL SYS	40-44
INTERLOCKING BOLT SYS	45-47
PROSTHESIS SYSTEM	48
SPINE SYSTEM	49-57



SMALL COMPRESSION LOCKING PLATE 3.5MM

Catalog no.	011.155.03 to 011.155.16
Screw ø	3.5mm
Material	Stainless Steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP SMALL T PLATE (OBLIQUE) 3.5MM - LEFT AND RIGHT

Catalog no.	011.151.02 to 011.151.10
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP SMALL T PLATE(RIGHT ANGLE) 3.5MM

Catalog no.	011.150.02 to 011.150.10
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)

LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 1)

Catalog no.	011.156.04 to 011.156.12
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,6,8,10,12
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 2)

Catalog no.	011.1910.03 to 011.1910.10
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 3)

Catalog no.	011.1911.03 to 011.1911.08
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, Titanium
Holes	3,4,5,6,7,8
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)





LCP CLOVERLEAF PLATE 3.5MM

Catalog no.	011.154.03 to 011.154.03
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP RECONSTRUCTION PLATE 3.5MM ROUND HOLE

Catalog no.	011.153.03 to 011.153.18
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)




LCP RECONSTRUCTION PLATE 5.0MM (ROUND HOLE)

Catalog no.	011.1954.03 to 011.1954.18
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm)

LCP RECONSTRUCTION PLATE 3.5MM (DYNAMIC HOLE)

Catalog no.	011.1912.03 to 011.1912.16
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP RECONSTRUCTION PLATE 5.0MM (DYNAMIC HOLE)

Catalog no.	011.1913.03 to 011.1913.16
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm)



DISTAL TIBIA LOCKING PLATE LEFT & RIGHT

Catalog no.	011.157.04 to 011.157.12
Screw ø	5.0mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4, 6, 8, 10, 12,
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)





MEDIAL DISTAL TIBIA LOCKING PLATE 3.5MM LEFT AND RIGHT

Catalog no.	011.1909.05 to 011.1909.13
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	5, 7, 9, 11, 13
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



DISTAL TIBIA ANATOMICAL LOCKING PLATE - LEFT & RIGHT

Catalog no.	011.1001.04 to 011.1001.14
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4, 6, 8, 10, 12, 14
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



NARROW LOCKING COMPRESSION PLATE 5.0MM

Catalog no.	011.166.04 to 011.166.18
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)

BROAD LOCKING COMPRESSION PLATE 4.8MM

Catalog no.	011.167.04 to 011.167.20
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
Shaft Part	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



L BUTTRESS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.148.04 to 011.148.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12
Proximal Screw	
Screw Length	
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



T PLATE LOCKING 5.0MM

Catalog no.	011.162.04 to 011.162.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12
Proximal Screw	
Screw Length	
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)





T BUTTRESS LOCKING PLATE 5.0MM

Catalog no.	011.1262.04 to 011.1262.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12
Proximal Screw	
Screw Length	
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



DISTAL HUMERUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1925.03 to 011.1925.15
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,15
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)



DISTAL HUMERUS WITH LATERAL SUPPORT LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1907.03 to 011.1907.15
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,15
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)

EXTRA-ARTICULAR DISTAL HUMERUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1941.04 to 011.1941.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,6,8,10,12,14
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)



LATERAL TIBIAL LOCKING PLATE 5.0MM/6.5MM – LEFT & RIGHT

Catalog no.	011.161.03 to 011.161.12
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Proximal Screw	Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locking Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



LATERAL TIBIAL LOCKING PLATE 5.0MM – LEFT & RIGHT

Catalog no.	011.159.03 to 011.159.13
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 5.0mm
Screw Length	30mm to 80mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locking Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)





DISTAL FEMUR LOCKING PLATE 5.0MM/6.5MM (LEFT AND RIGHT)

Catalog no.	011.163.03 to 011.163.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



DISTAL FEMUR LOCKING PLATE LEFT & RIGHT

Catalog no.	011.164.03 to 011.164.13
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 5.0mm
Screw Length	30mm to 80mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locking Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



CONDYLAR BUTTERS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.176.03 to 011.176.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locking Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)

DHS PLATE LOCKING LONG BARRAL 130° /135°

Catalog no.	011.165.02 to 011.165.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)



DHS PLATE LOCKING SHORT BARRAL 130° /135°

Catalog no.	011.152.02 to 011.152.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)



D.C.S PLATE LOCKING 90°

Catalog no.	011.1937.03 to 011.1937.16
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Locking Screw 5.0mm	10mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)
Cancellous Locking Screw 6.5mm	30mm to 95mm (Diff. 5.0mm)





DISTAL FEMUR LOCKING PLATE 5.0MM/6.5MM (LEFT AND RIGHT)

Catalog no.	011.1928.04 to 011.1928.08
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	04, 06, 08
Locking Screw	
3.5mm	12mm to 40mm (Diff. 2.0mm)



PROXIMAL FEMUR LOCKING PLATE

Catalog no.	011.1938.05 to 011.1938.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	05, 07, 09, 11, 13
Locking Screw	
5.0mm	10mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)
Cancellous	
Screw 6.5mm	30mm to 100mm (Diff. 5.0mm)



PROXIMAL TIBIA LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1986.04 to 011.1986.14
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	04, 06, 08, 10, 12, 14
Locking Screw	
3.5mm	12mm to 40mm (Diff. 2.0mm)


LCP TOMIFIX MID HIGH T PLATE

Catalog no.	011.1252.03 & 011.1252.04
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	03, 04
Locking Screw	
5.0mm	18mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)



CALCANEUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1940.70mm & 011.1940.75mm
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	70mm, 75mm
Locking Screw	
3.5mm	12mm to 40mm (Diff. 2.0mm)



DISTAL TIBIA WITH TAB LOCKING PLATE - LEFT & RIGHT

Catalog no.	011.1967.05 to 011.1967.13
Screw ø	4.0 mm & 3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	05, 07, 09, 11, 13
Locking Screw	
3.5mm	10mm to 40mm (Diff. 2.0mm)
Locking Screw	
4.0mm	20mm to 70mm (Diff. 5.0mm)



2.4MM DISTAL RADIUS LOCKING PLATE STRAIGHT

Catalog no.	011.193.05 to 011.193.13
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	5, 7, 9, 11, 13
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



2.4MM L DISTAL RADIUS LOCKING PLATE (HEAD 2HOLE) – LEFT & RIGHT

Catalog no.	011.191.3 & 011.191. 4
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2+3, 2+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



2.4MM L DISTAL RADIUS LOCKING PLATE (HEAD 3HOLE) – LEFT & RIGHT

Catalog no.	011.190. 3 & 011.190. 4
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3+3, 3+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



2.4MM L DISTAL RADIUS OBLIQUE LOCKING PLATE – LEFT & RIGHT

Catalog no.	011.194. 3 & 011.194. 4
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3+3, 3+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)




2.4MM T DISTAL RADIUS LOCKING PLATE (HEAD 3HOLE)

Catalog no.	011.192. 3 & to 011.192. 4
Screw ø	2.4mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3+3, 3+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



2.4MM DISTAL RADIUS LOCKING PLATE (HEAD 5HOLE) – LEFT & RIGHT

Catalog no.	011.197. 3 & to 011.197. 5
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	5+3, 5+5
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)





2.4 MM BUTTRESS LOCKING PLATE (HEAD 5HOLE) – LEFT & RIGHT

Catalog no.	011.198
Screw ø	2.4mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	5+3
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



TWO-COLUMN DISTAL RADIUS PLATE LEFT & RIGHT

Catalog no.	011.1946.02 to 011.1946.04
Screw ø	2.4mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	02, 03, 04
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)

LCP SCREW DIA 2.4 MM

Catalog no	011.286.06 to 011.286.30
Thread Diameter	2.4mm
Use drills Diameter	2.0mm
Drive	hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)




LCP SCREW DIA 2.7 MM

Catalog no	011.204.06 to 011.204.30
Thread Diameter	2.7mm
Use drills Diameter	2.0mm
Drive	hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)



LCP SCREW DIA 3.5 MM

Catalog no	011.268.10 to 011.268.50
Thread Diameter	
Use drills Diameter	2.7mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	3.5mm
Length	10mm to 50mm (Diff. 2.0mm)





LCP SCREW DIA 4.0 MM

Catalog no	011.270.25 to 011.270.70
Thread Diameter	4.0mm
Use drills Diameter	3.2mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	4.0mm
Length	25mm to 70mm (Diff. 2.0mm)



LCP SCREW DIA 5.0 MM

Catalog no	011.269.16 to 011.269.95
Thread Diameter	5.0mm
Use drills Diameter	4.0mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	5.0mm
Length	16mm to 50mm (Diff. 2.0mm) 55mm to 95mm (Diff. 5.0mm)



LCP SCREW DIA 6.5 MM

Catalog no	011.271.35 to 011.271.95
Thread Diameter	6.5mm
Use drills Diameter	4.5mm
Drive	hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	6.5mm
Length	35mm to 95mm (Diff. 5.0mm)

RECONSTRUCTION PLATE FOR 2.7MM SCREW ROUND HOLE

Catalog no.	011.300.04 to 011.300.14
Screw ø	2.7mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 2.7mm
Screw Length	06mm to 30mm (Diff. 2.0mm)



DYNAMIC COMPRESSION PLATE FOR 2.7 MM SCREW

Catalog no.	011.320.04 to 011.320.14
Screw ø	2.7mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 2.7mm
Screw Length	06mm to 30mm (Diff. 2.0mm)



ONE THIRD TUBULAR PLATE FOR 3.5 MM SCREW

Catalog no.	011.321.04 to 011.321.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)





RECONSTRUCTION PLATE FOR 3.5 MM SCREW ROUND HOLE

Catalog no.	011.322.04 to 011.322.18
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



RECONSTRUCTION PLATE FOR 3.5 MM SCREW DC HOLE

Catalog no.	011.3221.04 to 011.3221.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



RECONSTRUCTION CURVED PLATE FOR 3.5MM SCREW

Catalog no.	011.301.04 to 011.301.16
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,08,09,10,11,12,13,14,15,16
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)

SMALL DCP PLATE FOR 3.5 MM SCREW

Catalog no.	011.325.04 to 011.325.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



LC DCP PLATE SMALL FOR 3.5 MM SCREW

Catalog no.	011.333.04 to 011.333.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



SMALL T PLATE FOR 3.5 MM SCREW

Catalog no.	011.326.03 to 011.326.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	03,04,05,06,07,08,09,10
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)





SMALL T OBLIQUE PLATE FOR 3.5 MM SCREW

Catalog no.	011.327.03 to 011.327.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	03,04,05,06,07,08,09,10
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



ACRO PLATE FOR 3.5 MM SCREW LEFT & RIGHT

Catalog no.	011.3931.04 to 011.3931.06
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



RECONSTRUCTION PLATE FOR 4.5 MM SCREW ROUND HOLE

Catalog no.	011.323.04 to 011.323.18
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)

RECONSTRUCTION PLATE FOR 4.5 MM SCREW DC HOLE

Catalog no.	011.324.04 to 011.324.14
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14,
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)




SEMI TUBULAR PLATE FOR 4.5 MM SCREW

Catalog no.	011.3225.04 to 011.3225.14
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14,
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



NARROW DCP PLATE FOR 4.5 MM SCREW

Catalog no.	011.337.04 to 011.337.14
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)





LC DCP PLATE NARROW FOR 4.5 MM SCREW

Catalog no.	011.334.04 to 011.334.14
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



BROAD DCP PLATE FOR 4.5 MM SCREW

Catalog no.	011.338.04 to 011.338.18
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



LC DCP PLATE BROAD FOR 4.5 MM SCREW

Catalog no.	011.358.06 to 011.358.14
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)

CLOVERLEAF PLATE

Catalog no.	011.328.03 to 011.328.12
Screw ø	
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Screw	
Screw Length	



T PLATE FOR 4.5 MM SCREW

Catalog no.	011.329.03 to 011.329.12
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



T BUTTRESS PLATE FOR 4.5 MM SCREW

Catalog no.	011.330.03 to 011.330.12
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)





L BUTTRESS PLATE FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.332.03 to 011.332.12
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



LATERAL TIBIAL PLATE (HOCKEY) FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.335.04 to 011.335.14
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



CONDYLAR BUTTRESS PLATE FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.336.04 to 011.336.14
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)

DHS PLATE (DCP HOLES) SHORT BARREL

Catalog no.	011.339.04 to 011.339.14
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Degree	120°, 125°, 130°, 135°, 140°, 145°
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



DHS PLATE (ROUND HOLES) SHORT BARREL

Catalog no.	011.339.04R to 011.339.18R
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Degree	120°, 125°, 130°, 135°, 140°, 145°
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



DHS PLATE (DCP HOLES) LONG BARREL

Catalog no.	011.340.04 to 011.340.14
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Degree	120°, 125°, 130°, 135°, 140°, 145°
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)





DHS PLATE (ROUND HOLES) LONG BARREL

Catalog no.	011.340.02R to 011.340.18R
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Degree	120°, 125°, 130°, 135°, 140°, 145°
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



DHS PLATE (ROUND HOLES) PEDIATRICS

Catalog no.	011.331.04 to 011.331.10
Screw ø	
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10
Degree	125°, 130°, 135°
Proximal Screw	Lag Screw
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



95° DCS PLATE (DCP HOLES) TOTAL HOLE

Catalog no.	011.341.06 to 011.341.20
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)

DCS PLATE (ROUND HOLES) 95°

Catalog no.	011.342.06 to 011.342.20
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)



Y RECONSTRUCTION PLATE FOR 3.5mm SCREW

Catalog no.	011.302.01 to 011.302.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	1,2,3,4,5,6,7,8,9,10
Distal Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	10mm to 30mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	35mm to 70mm (Diff. 5.0mm)



CALCANEAL PLATE 3.5MM

Catalog no.	011.310
Screw ø	4.0mm & 3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	Small , Medium , Large
Distal Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	10mm to 30mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm





COBRA HEAD PLATE FOR 4.5mm SCREW LEFT & RIGHT

Catalog no.	011.311.04 to 011.311.14
Screw ø	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)



TROCENTRIC PLATE FOR DHS

Catalog no.	011.3948
Screw ø	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4 Hole
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)



EIGHT PLATE

Catalog no.	011.366
Screw ø	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4 Hole
Distal Screw	Cortical Screw 4.5mm

CORTICAL SCREW (SELF TAPPING)

1.5MM

Catalog no.	011.49.06 to 011.49.30
Material	Stainless steel 316L,316LVM, Titanium
Size	1.5mm
Length	06mm to 30mm (Diff. 2.0mm)



CORTICAL SCREW

2.0MM

Catalog no.	011.49.06 to 011.49.30
Material	Stainless steel 316L,316LVM, Titanium
Size	1.5mm
Length	06mm to 30mm (Diff. 2.0mm)



CORTICAL SCREW

2.4MM

Catalog no.	011.485.06 to 011.485.30
Material	Stainless steel 316L,316LVM, Titanium
Size	2.4mm
Length	06mm to 30mm (Diff. 2.0mm)





CORTICAL SCREW 2.7MM

Catalog no.	011.401.06 to 011.401.30
Material	Stainless steel 316L,316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)



CORTICAL SCREW (14 TPI) 3.5MM

Catalog no.	011.402.10 to 011.402.30
Material	Stainless steel 316L,316LVM, Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



CORTICAL SCREW (14 TPI) SELF TAPPING 3.5MM

Catalog no.	011.01.402.10 to 011.01.402.60
Material	Stainless steel 316L,316LVM, Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)

CORTICAL SCREW 3.5MM (20 TPI)

Catalog no.	011.403.10 to 011.403.10.60
Material	Stainless steel 316L,316LVM Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



CORTICAL SCREW (20 TPI) SELF TAPPING 3.5MM

Catalog no.	011.01.403.10 to 011.01.403.60
Material	Stainless steel 316L,316LVM, Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



CANCELLOUS SCREW (SHORT THREAD) 4.0MM

Catalog no.	011.404.10 to 011.404.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm) 35mm to 70mm(Diff. 5.0mm)





CANCELLOUS SCREW (FULL THREAD) 4.0MM

Catalog no.	011.405.10 to 011.405.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm) 35mm to 70mm(Diff. 5.0mm)



CANNULATED SCREW (SHORT THREAD) 4.0MM

Catalog no.	011.406.10 to 011.406.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm) 35mm to 70mm(Diff. 5.0mm)



CANNULATED SCREW (FULL THREAD) 4.0MM

Catalog no.	011.407.10 to 011.407.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm) 35mm to 70mm(Diff. 5.0mm)

CORTICAL SCREW 4.5MM DIA

Catalog no.	011.409.16 to 011.409.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	16mm to 70mm(Diff. 2.0mm)



CORTICAL SCREW (SELF TAPPING) 4.5MM

Catalog no.	011.01.409.12 to 011.01.409.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	16mm to 70mm(Diff. 2.0mm)



MALLEOLAR SCREW 4.5MM

Catalog no.	011.410.20 to 011.410.100
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	20mm to 100mm(Diff. 5.0mm)



CANCELLOUS SCREW (16TH) 6.5MM

Catalog no.	011.411.20 to 011.411.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)



CANCELLOUS SCREW (32TH) 6.5MM

Catalog no.	011.412.20 to 011.412.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)



CANCELLOUS SCREW (FULL TH) 6.5MM

Catalog no.	011.413.25 to 011.413. 110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	25mm to 110mm(Diff. 5.0mm)



CANNULATED SCREW (16TH) 6.5MM

Catalog no.	011.414.20 to 011.414.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)



CANNULATED SCREW (32TH) 6.5MM

Catalog no.	011.415.25 to 011.415.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	25mm to 110mm(Diff. 5.0mm)



CANNULATED SCREW (FULL THREAD) 6.5MM

Catalog no.	011.417.20 to 011.417.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 100mm(Diff. 5.0mm)





HARBAT SCREW

Catalog no.	011.416.25 to 011.416.30
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	25mm to 30mm(Diff. 2.0mm)



DHS LAG SCREW

Catalog no.	011.418.50 to 011.418.110
Material	Stainless steel 316L,316LVM, Titanium
Size	12.5mm
Length	50mm to 110mm(Diff. 5.0mm)



DHS LAG SCREW (HEX)

Catalog no.	011.01.418.50 to 011.01.418.110
Material	Stainless steel 316L,316LVM, Titanium
Size	12.5mm
Length	50mm to 110mm(Diff. 5.0mm)

DHS LAG SCREW PEDIATRICS

Catalog no.	011.463.40 to 011.463.105
Material	Stainless steel 316L,316LVM, Titanium
Size	
Length	40mm to 105mm(Diff. 5.0mm)



WASHER FOR CANCELLOUS SCREW

Catalog no.	011.4216
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm, 7.0mm, 13mm
Length	Standard



K NAIL FOR FEMUR

Catalog no.	011.545.20 to 011.545.48
Locking Bolt ø	
Material	Stainless steel 316L,316LVM ,Titanium
Nail Diameter	07mm, 08mm, 09mm, 10mm, 11mm, 12mm
Nail Length	20cm to 46cm (Diff. 2.0cm)



V NAIL FOR TIBIA

Catalog no.	011.546.20 to 011.546.40
Locking Bolt ø	
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	07mm, 08mm, 09mm, 10mm, 11mm,
Nail Length	20cm to 40cm (Diff. 2.0cm)



SUPRA CONDYLAR NAIL

Catalog no.	011.548.15 to 011.548.30
Locking Bolt ø	
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	15cm to30cm(diff.5.0CM)
Proximal & Distal Bolt Bolt Length	



RUSH NAIL

Catalog no.	011.534.18 to 011.534.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	1.5mm, 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm, 6mm
Nail Length	18cm to 30cm (diff. 2.0cm)



SQUARE NAIL FOR ULNA

Catalog no.	011.547.2.0.18 to 011.547.2.0.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)



SQUARE NAIL FOR RADIUS

Catalog no.	011.547.2.0.18 to 011.547.2.0.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)





HUMERUS INTERLOCKING NAIL

Catalog no.	011.549.18 to 011.549.30
Locking Bolt ø	3.1mm
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	6.0mm, 7.0mm, 8.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)
Proximal Bolt	Locking Bolt 3.1mm
Bolt Length	10mm to 40mm(diff. 2.0mm)
Distal Bolt	Locking Bolt 3.1mm
Bolt Length	10mm to 40mm(diff. 2.0mm)



INTER LOCKING TIBIA NAIL (SOLID)

Catalog no.	011.01.550.26 to 011.01.550.40
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM ,Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm
Nail Length	26.0cm to 40.0cm(diff. 2.0cm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)



INTER LOCKING TIBIA NAIL

Catalog no.	011.550.26 to 011.550.40
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM ,Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm
Nail Length	26.0cm to 40.0cm(diff. 2.0cm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)

INTER LOCKING FEMUR NAIL

Catalog no.	011.551.34. to 011.551.44
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm,12.0mm
Nail Length	34.0cm to 44.0cm (Diff. 2.0cm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)



TFN NAIL SHORT

Catalog no.	011.5935
Locking Bolt ø	8.0mm,6.4mm & 4.9mm
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	9.0mm, 10.0mm, 11.0mm, 12.0mm
Nail Length	18cm
Degree	130° /135°
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm(diff. 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm(diff. 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)



PFN NAIL SHORT

Catalog no.	011.57
Locking Bolt ø	5.0mm, 3.7mm & 3.0mm
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	9.0mm, 10.0mm, 11.0mm, 12.0mm
Nail Length	25cm
Degree	130° /135°
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm(diff. 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm(diff. 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)






ELASTIC NAIL (FLEXIBLE NAIL)

Catalog no.	011.5938
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm
Nail Length	44.0cm Standard



SIRUS NAIL (LEFT & RIGHT)

Catalog no.	011.559.34 to 011.559.44
Locking Bolt ø	6.5mm & 4.9mm
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	9.0mm,10.0mm,11.0mm,12.0mm
Nail Length	34cm to 44cm (diff. 2.0cm)
Proximal Bolt	Cannulated / Solid 6.5mm
Bolt Length	30mm to 110m (diff. 5.0mm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm) 65mm to 80mm(diff. 5.0mm)
Standard	Left & Right



ENDER'S NAIL

Catalog no.	011.553.18 to 011.553.46
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	3.0mm, 3.5mm, 4.0mm, 4.5mm,
Nail Length	18cm to 46cm (Diff. 1.0cm)

INTER LOCKING BOLT 3.1MM

Catalog no.	011.633.18 to 011.633.50
Thread Diameter	3.1mm
Material	Stainless steel 316L,316LVM ,Titanium
Size	3.1mm
Length	18mm to 50mm (Diff. 2.0mm)



INTER LOCKING BOLT 4.5MM

Catalogue No.	011.01.652.18 to 011.01.652. 90
Thread Diameter	4.5mm
Core Diameter	
Use drills Diameter	3.2mm
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	18mm to 90mm (Diff. 2.0mm)



INTER LOCKING BOLT 3.9MM

Catalog no.	011.653.18 to 011.653.80
Thread Diameter	3.9mm
Core Diameter	2.7mm
Use drills Diameter	3.2mm
For passing hole	4.0mm
Material	Stainless steel 316L,316LVM, Titanium
Size	3.9mm
Length	18mm to 50mm (Diff. 2.0mm) 55mm to 80mm (Diff. 5.0mm)





INTER LOCKING BOLT 4.9MM

Catalog no.	011.652.18 to 011.652.90
Thread Diameter	4.9mm
Core Diameter	3.8mm
Use drills Diameter	4.0mm
For passing hole	5.1mm
Material	Stainless steel 316L,316LVM, Titanium
Size	4.9mm
Length	18mm to 60mm (Diff. 2.0mm) 65mm to 90mm (Diff. 5.0mm)

PFN BOLT CANNULATED 6.4MM

Catalog no.	011.65.30 to 011.65.120
Thread Diameter	
Core Diameter	
Use drills Diameter	
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	6.4mm
Length	30mm to 120mm (Diff. 5.0mm)



PFN BOLT CANCELLOUS 6.4MM

Catalog no.	011.01.65.30 to 011.01.65.120
Thread Diameter	
Core Diameter	
Use drills Diameter	
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	6.4mm
Length	30mm to 120mm (Diff. 5.0mm)



PFN BOLT CANNULATED 8.0MM

Catalog no.	011.66.45 to 120
Thread Diameter	
Core Diameter	
Use drills Diameter	
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	8.0mm
Length	45mm to 120mm (Diff. 5.0mm)



PFN BOLT CANNULATED 8.0MM (HEX)

Catalog no.	011.01.66.45 to 011.01.66.120
Thread Diameter	
Core Diameter	
Use drills Diameter	
For passing hole	
Material	Stainless steel 316L,316LVM ,Titanium
Size	8.0mm
Length	45mm to 120mm (Diff. 5.0mm)





THOMPSON HIP PROSTHESIS

Catlog no.	011.754.37 to 011.754.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Standard



AUSTIN MOORE HIP PROSTHESIS

Catlog no.	011.755.37 to 011.755.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Standard



SHOULDER PROSTHESIS

Catlog no.	011.751.35 to 011.751.51
Material	Stainless steel 316L, 316LVM, Titanium
Size	35, 37, 39, 41, 43, 45, 47, 49, 51
Length	Standard

MONO AXIAL SCREW

Catalog no.	011.826.20 to 011.826.50
Material	Titanium
Dia	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



POLY AXIAL SCREW

Catalog no.	011.828.20 to 011.828.50
Material	Titanium
Dia	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



MONO AXIAL SCREW (D/T)

Catalog no.	011.833.20 to 011.833.50
Material	Titanium
Dia	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)





POLY AXIAL SCREW (D/T)

Catalog no.	011.834.20 to 011.834.50
Material	Titanium
Dia	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



MONO AXIAL REDUCTION SCREW

Catlog no.	011.827.20 to 011.827.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



POLY AXIAL REDUCTION SCREW

Catlog no.	011.829.20 to 011.829.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)

MONO AXIAL REDUCTION (D/T)

Catlog no.	011.813.20 to 011.813.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



POLY AXIAL REDUCTION (D/T)

Catlog no.	011.895.20 to 011.895.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



3.5 MM DIA - MONOAXIAL SCREW

Catalog no.	011.830.12 to 011.830.40
Material	Titanium
Dia	3.5mm, 4.2mm, 5.0mm, 6.0mm
Length	12mm to 20mm(Diff. 2.0mm) 25mm to 40mm(Diff. 5.0mm)





3.5MM DIA – POLYAXIAL SCREW

Catalog no.	011.831.12 to 011.831.40
Material	Titanium
Dia	3.5mm, 4.2mm
Length	12mm to 20mm(Diff. 2.0mm) 25mm to 40mm(Diff. 5.0mm)



LAMINAR HOOK ANGLE BLADE

Catalog no.	011.835.01
Material	Titanium
Dia	Standard
Length	Standard



LAMINAR HOOK NARROW BLADE

Catalog no.	011.831.12 to 011.831.40
Material	Titanium
Dia	3.5mm, 4.2mm
Length	12mm to 20mm(Diff. 2.0mm) 25mm to 40mm(Diff. 5.0mm)

PEDICAL HOOK

Catlog no.	011.835.03
Material	Titanium
Size	Standard
Length	Standard



LEMNAR HOOK WIDE BLADE

Catlog no.	011.835.04
Material	Titanium
Size	Standard
Length	Standard



STAPLE

Catlog no.	011.829.01
Material	Titanium
Size	Standard
Length	Standard





CONNECTING ROD

Catalog no.	011.832
Material	Titanium
Dia	3.0mm, 3.5mm, 5.0mm, 5.5mm
Length	50mm to 450mm



3.5mm Dia - POSTERIOR OCCIPITAL PLATE WITH ROD

Catalog no.	011.854.3.5.100 to 011.854.3.5.300
Material	Titanium
Dia	3.5mm
Length	100mm, 200mm, 250mm, 300mm



POSTERIOR OCCIPITAL SCREW

Catalog no.	011.8133.06 to 011.8133.20
Material	Titanium
Dia	3.5mm
Length	06mm to 20mm (Diff. 2.0mm)

ANTERIOR CERVICAL CAGE (SQUARE TYPE)

Catlog no.	011.844.4 to 8.14 to 11.844.4 to 8.16
Material	Titanium
Size	4.0mm, 5.0mm, 06mm, 7.0mm, 8.0mm
Length	12mmx14mm, 12mmx16mm



EXPANDABLE CAGE

Catlog no.	011.846.20 to 011.846.75
Material	Titanium
Size	10mm, 12mm, 14mm, 16mm, 18mm, 20mm
Length	20mm to 75mm (Diff. 5.0mm)



ECO CAGE

Catlog no.	011.845.15 to 011.845.100
Material	Titanium
Size	10mm, 12mm, 14mm, 16mm, 18mm, 20mm
Length	15mm to 100mm (Diff. 5.0mm)





PLIF CAGE

Catalog no.	011.847.20 to 011.847.30
Material	Titanium
Size	7.0mm, 8.0mm, 9.0mm, 10.0mm, 11.0mm
Length	20mm, 25mm, 30mm



TLIF CAGE (BANANA TYPE)

Catalog no.	011.856.7 to 011.856.11
Material	Titanium
Dia	7.0mm, 8.0mm, 9.0mm, 10.0mm, 11.0mm
Length	4o, 8o, 11o



PLIF CAGE (BULLET TYPE)

Catalog no.	011.857.20 to 011.857.30
Material	Titanium
Dia	8.0mm, 10.0mm, 12.0mm
Length	20mm, 25mm, 30mm

ANTERIOR CERVICAL PLATE (Z TYPE)

Catalog no.	011.842.25 to 011.842.70
Material	Titanium
Dia	Standard
Length	25.0mm,27.5mm,30.0mm,32.5mm, 35.0mm,37.5mm,40.0mm,42.5mm, 45.0mm,47.5mm,50.0mm,55.0mm, 60.0mm,65.0mm,70.0mm



BONE SCREW FOR ACP 4mm TITANIUM

Catalog no.	011.843.10 to 011.843.24
Material	Titanium
Dia	4.0mm
Length	10mm to 24mm (Diff. 2.0mm)



OCCIPITAL LOOPS

Catlog no.	011.855
Material	Titanium
Size	4.0mm
Length	



CONTACT US

A great vision without great people is irrelevant. Let's cooperate.

📍 Jouriyan Kalan, Sialkot 51310 Pakistan.

☎ +92-300-614 11 83

📞 (052) 355-87-36

✉ info@lasanisurgicalind.com

🌐 www.lasanisurgicalind.com